

## CRUTCHES HOLDER AND CLAMP FOR MOBILITY AID



### **Description:**

Trying to carry crutches, gear, or everyday items at the same time is not always practical. This adjustable clamp gives you a simple way to secure your crutches or attach gear to your wheelchair, walking frame, bike, or frame so you can move more freely without juggling everything at once. The dual clamp design connects two items securely, helping keep things in place while still being easy to adjust or remove when you need to. It is a practical add-on that makes everyday movement, transport, and transitions between supports a lot more manageable.

### **Great for:**

-  Attaching crutches to a wheelchair or walking frame
-  Carrying crutches while using another support
-  Securing items to bike or scooter handlebars
-  Keeping two crutches clipped together
-  Freeing up your hands while moving
-  Mounting a torch or light to your chair or frame
-  Holding a lightweight umbrella in place
-  Holding canes or similar mobility aids
-  Transporting items without juggling them

### **Features:**

-  Dual clamp holds two items securely
-  Adjustable to fit different tube sizes
-  Quick to attach and remove
-  Strong plastic and steel build

**Details:**

**Material:** Plastic and Steel

**Size:** 10cm (L) x 5.6cm (W) x 2cm (H) (centre twist base is 3.2cm wide)

*Fits item with a diameter between 1.8 - 3.8cms.*

**Weight:** Approx 80 grams

**Package Includes:** 1x Crutches Clamp

**Age Guidance:** Suitable for teens and adults; younger users may need assistance when attaching or adjusting. *Not intended as a toy.*

**How to Use:**

*Attach:*

 Clamp one side onto your wheelchair, frame, or base

 Tighten so it stays in place

*Clip in:*

 Open the other side to add your crutch or item

 Adjust so it sits comfortably

*Use:*

 Clip items in and out as needed

 Check it's secure before moving

**Care Instructions:**

 Wipe clean with a damp cloth

 Dry fully after cleaning

 Avoid harsh chemicals or soaking