



CAR STANDING TRANSFER AID

INSTRUCTION MANUAL



Getting in and out of a vehicle can take more movement, balance, and effort than people often realise, especially with taller vehicles, lower seats, long days, busy outings, or when a little extra support feels helpful.

This car standing transfer aid creates an additional grab handle inside the vehicle to assist with entering, exiting, repositioning, or standing up from the seat. The adjustable strap loops securely around the door frame, creating a sturdy support point that is easy to grab when needed.

Lightweight and easy to keep in the car, it can be helpful for everyday travel, accessibility support, appointments, road trips, recovery periods, or simply making vehicle transfers feel more comfortable and supported.

Great for:

- Getting in and out of vehicles more comfortably
- Hip, knee, or leg injuries and recovery
- Balance or stability support during transfers
- Extra support when standing or repositioning
- Accessibility and mobility support
- Wheelchair users and transfer assistance
- Surgery recovery or rehabilitation periods
- Spinal cord injury and neurological support needs
- Low-energy or fatigue-heavy days
- SUVs, vans, utes, and taller vehicles
- Everyday errands, outings, and appointments
- Added confidence during car transfers
- Independent or assisted vehicle access
- Keeping in the car for everyday convenience
- Long travel days and repeated car transfer

Features:

- Creates an extra grab handle inside the vehicle
- Adjustable strap attachment design
- Supports standing, repositioning, and transfers
- Compatible with most vehicle door frames
- Soft non-slip grip handle
- Lightweight and portable
- Easy to store in the car between uses
- Supports up to 100kg



CAR STANDING TRANSFER AID

INSTRUCTION MANUAL

Image Description: A black vehicle transfer support handle attached to the upper door frame inside a red car. The adjustable strap creates a hanging grab handle to assist with entering, exiting, standing, or repositioning in the vehicle. Additional images show the handle attached in different positions and being used during transfers.

Details:

Material: ABS plastic, canvas strapping

Size: 41cm (L) x 12cm (W) x 1cm (D)

Weight: 120 grams

Package Contains: 1x Pull Assist Handle Transfer Aid

Age Guidance: Suitable for teens and adults. Younger users may use with support where appropriate.

How to Use:

1. Loop the strap securely around the vehicle door frame
2. Adjust strap length if needed
3. Use the handle for support when entering or exiting the vehicle
4. Use during standing, lowering, or repositioning as needed
5. Can be left installed, stored in the vehicle, or taken with you between cars

Care Instructions:

- Wipe clean with a damp cloth
- Allow to dry fully before storing
- Store in a clean, dry place
- Check strap and handle regularly for wear or damage

Please Note:

- Maximum safe supported weight is 100kg
- Ensure the strap is securely fitted before use
- Do not use if the strap or handle shows signs of damage